

# CooperSurgical

## HBA®

### Sperm-Hyaluronan Binding Assay

CE  
RxOnly

#### PROPRIETARY AND ESTABLISHED NAMES:

HBA® Sperm-Hyaluronan Binding Assay

#### INTENDED USE:

The HBA® Sperm-Hyaluronan Binding Assay is an in vitro diagnostic indicated for use:

1. As a component of the standard analysis of semen in the diagnosis of suspected male infertility.
2. As a component of analyses of either raw or processed semen for determining the proper course of IVF treatment of infertility.

#### SUMMARY:

In natural fertilization mature sperm bind to hyaluronan, the main component of the *cumulus oophorus* matrix; immature sperm do not bind (3, 9, 13, 14, 19, 21).

Mature sperm also bind to hyaluronan chemically attached to a support (9, 10) such as the hyaluronan-coated glass slides that make up the HBA® Binding Assay. Viewed in the microscope, bound sperm are differentiated from unbound sperm by their beating tails with heads that make no progressive movement.

Hyaluronan binding indicates the successful completion of spermiogenic events mediated by the HspA2 chaperone protein (1, 2, 8). In the absence of the chaperone, deficits in meiosis, DNA repair, hyaluronan binding and other steps in sperm maturation may occur.

Sperm unable to bind hyaluronan share many aspects of immaturity: they retain cytoplasm and histones (9), they show a higher frequency of aberrant morphology (18, 22, 23) and they have a lower genomic integrity than hyaluronan binders – hyaluronan binders have many fewer DNA single strand breaks than non-binders (24) and a 4-to-6 fold lower frequency of chromosomal aneuploidies (10, 11). Sperm motility is stimulated upon binding hyaluronan (6, 7, 14). Thus, the ability to bind hyaluronan designates a mature sperm and the proportion of mature sperm indicates the maturity of sperm in a semen sample.

#### THE HBA® SLIDE:

The HBA® slide has two identical hyaluronan-coated assay chambers. By counting the numbers of motile bound and unbound sperm in a common area, the proportion of binders can be calculated.

#### SPECIMEN COLLECTION AND PREPARATION:

Semen should be obtained by masturbation, preferably following 2-3 days of abstinence, and collected in a clean, dry container. Keep it at 20-28°C for 30 minutes to liquefy. Samples still viscous after 30 minutes can be diluted in an equal volume of a sperm dilution medium. Perform the assay on the semen specimen within 3 hours of collection.

#### WARNINGS AND PRECAUTIONS:

1. *Semen is a hazardous fluid and must be handled and disposed of using normal precautions for infectious, biohazardous materials. Wear gloves, eye protection and a laboratory coat. Disinfect all waste and dispose of it in a biohazard container.*

2. *Do not use HBA® slides for selecting sperm for ICSI or IVF! The slides are not sterile, they contain endotoxins and may be embryotoxic. Purchase of this kit does not allow sperm selection on HBA® slides.*
3. *Unused slides may be discarded in regular municipal trash containers. However, slides that have been exposed to body fluids must be discarded as bio-hazardous waste. Broken slides should be discarded as sharps.*

#### STORAGE:

HBA® slides should be stored in a clean, dry environment at 20-28°C.

#### PERFORMING THE ASSAY:

##### Materials Provided:

HBA® slides  
Cover slips  
Product Instructions for Use

##### Materials not Provided:

Microscope with 400x magnification  
Disposable, sterile capped test tubes  
Mechanical counting device  
Human tubal fluid (HTF) or other acceptable sperm diluent containing >0.50 mg/mL protein  
Disposable specimen handling gloves  
Timer/clock  
Pipette to deliver 2-10 µL, disposable pipette tips  
Microscope reticle delineating area

#### STEP-BY-STEP PROCEDURE:

1. Perform the assay at 20-30° C.
2. Load one of the two assay chambers with semen.  
  
Immediately before use, mix the semen and pipette a drop, 7-10 µL in volume, onto an area near the outside edge of the assay chamber, making certain it stays within the chamber.
3. Immediately install the cover slip: avoid entrapping air bubbles.

The cover slip provides a grid of 100 squares, 0.1 x 0.1 mm, within a viewing circle. Install the cover slip with the circle upward. Avoid air bubbles by contacting one edge of the cover slip to the edge of the chamber, forming an angle like two arms of a hinge (Figure 1); slowly lower the opposite edge of the cover slip, closing the hinge – when the cover slip is quite close to the slide, it will contact the sample drop (Figure 2); continue to lower the cover slip without releasing it; release the cover slip when the circle is completely covered with sample (Figure 3); finally center the circle over the chamber.

Figure 1

Contact cover slip to the slide:



Figure 2




Figure 3



4. Once the coverslip is placed, the incubation period begins. Do not delay placement of the coverslip. Uneven distribution of binding may occur, affecting results. Incubate the chamber for at least 10 and not more than 20 minutes. Count the unbound, motile sperm and the bound, motile sperm in the same number of grid squares (area).

In 10 minutes all the sperm will contact and bind to the immobilized hyaluronan layer. After 20 minutes weak motile sperm may begin to lose motility.

Bound motile sperm will cease progressive movement but retain active tail beating. Dead and non-motile sperm show no tail movement. Non-binding motile sperm swim about freely. Acceptable assay precision will be achieved when the sum of bound plus unbound motile sperm is greater than 200. With oligospermic samples, it may be necessary to count the sperm in several whole microscope fields, outside the grid area marked on the coverslip, in order to count 200 motile sperm.

It will not always be possible to find 200 motile sperm in a semen sample, particularly in oligospermic samples

5. Calculate the percent hyaluronan-binding sperm.

The percent of sperm binding to the hyaluronan layer is calculated as follows:

$$\% \text{ Bound} = 100 \times \frac{\text{Bound Motile Sperm}}{\text{Bound Motile Sperm} + \text{Unbound Motile Sperm}}$$

Example:

A count of 89 unbound, motile sperm was obtained by counting 100 grid squares. Immediately after, 27 bound, motile sperm were counted in the 100 grid squares.

$$23 \% \text{ Bound} = 100 \times \frac{27 \text{ Bound Motile Sperm}}{27 \text{ Bound Motile Sperm} + 89 \text{ Unbound Motile Sperm}}$$

## PRECISION OF THE ASSAY:

The precision of the assay is related to the number of sperm counted. The more sperm counted, the lower the variance. Preferably count at least 100 total motile sperm; results are acceptable when at least 30 total motile sperm are counted.

In a study of variability within and among three production lots of slides, single semen specimens were assayed in replicates of assay slides. With a normal specimen, counting approximately 196 sperm per replicate, the results (means (% bound)/CVs (%))

and numbers of replicates) were: (86.9, 85.5, 85.0/4, 5, 4/20, 10, 10). With an abnormal, oligospermic specimen, counting approximately 41 sperm, the results were: (59.0, 55.8, 58.4/18, 34, 17/20, 10, 10).

## QUALITY CONTROL:

Fresh mature semen can be used to demonstrate the sperm-hyaluronan binding reaction; however, fresh semen controls are not commercially available. Users are encouraged to establish their own sources of fresh control semen and to develop an associated history of normal and abnormal HBA® scores.

The sperm-hyaluronan binding reaction is not altered by freezing and thawing of sperm (16, 20, 26). Similarly, sperm processed either by swim up or by gradient preparation bind normally to hyaluronan.

## EXPECTED RESULTS:

The percent of hyaluronan-bound sperm ranges from essentially zero to essentially 100%. In a study of 157 semen samples from patients presenting for diagnosis of suspected infertility, the upper half of the samples showed an average HBA® score of  $93 \pm 2.6\%$  binding. The lower half had HBA® scores from 88% to zero binding, the average was  $70 \pm 18\%$  binding. Approximately 1% of these samples showed no detectable binding to hyaluronan.

Sperm preparation (swim-up, gradient) can alter the proportion of HA binders in a semen sample.

## LIMITATIONS AND TROUBLESHOOTING:

Use only fresh semen, less than three hours old.

The sample should permit at least 30 motile sperm to be counted. If less than 30 sperm occur in the grid, install a microscope reticle or count successive different microscope fields to permit inspection of an area large enough to evaluate at least 30 motile (bound plus unbound) sperm. Alternatively, concentrate the sample and re-assay the concentrate.

Always be sure that the sample has been applied to the coated side of the slide. The coated side is up when "HYDAK® HBA" can be clearly read. Application of sperm and cover slip to the uncoated side will give the appearance of no sperm binding or, if the cover slip has been pressed down hard, the sperm may appear immotile.

**Whenever no binding is observed, repeat the assay of the sample on a different slide to confirm the result.**

Occasionally defective cover slips occur with the grid etched on top of the slip instead of on the side in contact with the sample. When this occurs, repeat the assay with a new cover slip. Contact your distributor for replacements.

If an air bubble blocks counting within the counting grid, repeat the assay with a new chamber and cover.

## SIGNIFICANCE OF THE RESULTS:

A low HBA® score reflects a low proportion of mature sperm in the semen and therefore predicts infertility. Hyaluronan-binding sperm, in contrast, are competent in the interaction with the oocyte and are associated with high genomic integrity (11, 14, 21, 25), which improves the quality of the paternal contribution to the zygote. Hyaluronan-binding differentiates high and low functional integrity and fertilizing potential.

In prediction of fertility, it is not necessary that all sperm should be capable of binding hyaluronan, but the proportion of hyaluronan binders must reach an effective level in order to achieve a good likelihood of fertility. Binding above the effective level probably does not further improve fertility. Based on correlations with normal sperm morphology, the Hyaluronan Binding Assay level differentiating higher from lower expectation of fertility is estimated to be approximately 80%, Table 1

Table 1

HBA® Score Results Interpretation	
HBA® Score (% binding)	Interpretation
≥80 % binding	Normal maturity and physiological function
<80 % binding	Diminished maturity and physiological function

### PERFORMANCE CHARACTERISTICS:

A study comparing HBA® score results to normal sperm morphology (strict criteria) was carried out with semen samples from three IVF or andrology laboratories in three states (CT, CA and PA), with patients presenting for diagnosis of suspected infertility. At least 50 samples were obtained from each site and the HBA® assay performed on site according to the HBA® Instructions for Use. Morphology slides were prepared on site and sent to a central location where they were scored by a single reader.

The data were analyzed by a two-way contingency table using morphology (classified as subfertile if <5% normal morphology, otherwise fertile (4, 5, 15, 17)) and HBA® score (classified as subfertile if <80% binding, otherwise fertile). Significance was tested with the Pearson Chi-squared statistic using SPSS for Windows software. The analysis showed that HBA® score was significantly related to morphology for data from the combined data (all three sites), Table 2.

The combined data showed good specificity and positive predictive values, although sensitivity was only 40%. Therefore, at a cutoff of <80% binding, the HBA® assay detects less than half of semen samples with truly low morphology, but among those detected as low HBA® score there is a strong prediction (>80%) of poor morphology.

In a subsequent clinical trial with human outcomes, an HBA® Score of 65% was used as the key cutoff (24).

Table 2

Significance of HBA® Score to Morphology					
Study Site	N	Sensitivity	Specificity	Positive Predictive Value	p =
PA	50	43.2	76.9	84.2	0.198 (NS)
PA-r <sup>2</sup>	50	51.3	90.9	95.2	0.012
CA	52	54.2	78.6	68.4	0.015
CT	55	26.5	100	100	0.01
Combined	157	40.0	85.5	80.9	0.001















<sup>c</sup>Site data recalculated with cutoffs of <5% for morphology and <80% for HBA®.

### REFERENCES:

- Allen, J.W. et al. 1996. "Hsp70-2 is part of the synaptonemal complex in mouse and hamster spermatocytes". *Chromosoma* 104:414-21
- Dix, D.J. et al. 1996. "Targeted gene disruption of Hsp 70-2 results in failed meiosis, germ cell apoptosis and male infertility." *Proc. Natl. Acad. Sci. USA* 93:3264-8.
- Drahorad, J., J. Tesarik, D. Cechova and V. Vilim. 1991. Proteins and glycosaminoglycans in the intercellular matrix of the human cumulus-oophorous and their effect on conversion of proacrosin to acrosin. *J. Reprod. Fert.* 93:253- 262.
- Guizick, D.S., et al. 2001. Sperm Morphology, Motility and Concentration in Fertile and Infertile Men. *N. Engl. J. Med.* 345:1388-93.
- Gunalp, S. et al. 2001. A Study of Semen Parameters with Emphasis on Sperm Morphology in a Fertile Population: an Attempt to Develop Clinical Thresholds. *Hum. Reprod.* 16:110-14.
- Hamamah, S., et al. 1996. Identification of hyaluronic acid and chondroitin sulfates in human follicular fluid and their effects on human sperm motility and the outcome of in vitro fertilization. *Reprod. Nutr. Dev.* 36:43-52.
- Huszar, G., M. Willetts and M. Corrales. 1990. Hyaluronic acid (Sperm Select) improves retention of sperm motility and velocity in normospermic and oligospermic specimens. *Fertil. Steril.* 54:1127-34.
- Huszar, G. et al. 2000. "Putative creatine kinase M- isoform in human sperm is identified as the 70- kilodalton heat shock protein HspA2. *Biol. Reprod.* 63:925-32.
- Huszar, G., et al. 2003. Hyaluronic acid binding by human sperm indicates cellular maturity, viability and normal chromatin structure. *Fertil. Steril.* 79:1616-24.

- Huszar, G., et al. 2003. "Sperm testing by hyaluronic acid binding: andrological laboratory assessment and sperm selection for ICSI" in *Proc. of the Sero Symposium International Conference on Biotechnology of Human Reproduction: From Research to Clinical Application*, Torino, Italy, May 10-11, 2002, pps 149-157.
- Jakab, A. et al. 2005. "Intracytoplasmic sperm injection: a novel selection method for sperm with normal frequency of chromosomal aneuploidies." *Fertil. Steril.* 84(6): 1665-73.
- Johnston, J.B. 2000. A simple, nondestructive assay for bound hyaluronan. *J. Biomed. Mater. Res. (Appl. Biomater.)* 53:188-91.
- Kornovski, B.S., et al. 1994. The regulation of sperm motility by a novel hyaluronan receptor. *Fertil. Steril.* 61:935-40.
- Lin, Y., et al. 1994. A hyaluronidase activity of the sperm plasma membrane protein PH-20 enables sperm to penetrate the cumulus cell layer surrounding the egg. *J. Cell Biol.* 125:115-63.
- Menkveld, R., et al. 2001. Semen Parameters, Including WHO and Strict Criteria Morphology, in a Fertile and Subfertile Population: an Effort Towards Standardization of In-Vivo Thresholds. *Hum. Reprod.* 16:1165-71.
- Nijs, M. et al. 2009. "Influence of freeze-thawing on hyaluronic acid binding of human spermatozoa" *Reprod. Biomed. Online* 19(2):202-6.
- Ombelet, W. et al. 1997. Semen Parameters in a Fertile Versus Subfertile Population: a Need for Change in the Interpretation of Semen Testing. *Hum. Reprod.* 12:987- 93.
- Prinosilova, P. et al. 2009. "Selectivity of hyaluronic acid binding for spermatozoa with normal Tygerberg strict morphology" *Reprod. BioMed. Online* 18(2): 177-83.
- Ranganathan, S., A.K. Ganguly and K. Datta. 1994. Evidence for presence of hyaluronan binding protein on spermatozoa and its possible involvement in sperm function. *Mol. Reprod. Dev.* 38:69-76.
- Reed, M. et al. 2009. "Soy lecithin replaces egg yolk for cryopreservation of human sperm without adversely affecting postthaw motility, morphology, sperm DNA integrity or sperm binding to Hyaluronate. *Fertil. Steril.* 92(5):1787-90.
- Sabeur, K., et al. 1997. The PH-20 protein in human spermatozoa. *J. Androl.* 18:151-8.
- Worrlow, K.C. et al. 2003. "Statistically significant correlation between normal sperm morphology and sperm binding potential to an immobilized monolayer of hyaluronan using Hydak® HBA Sperm Binding Assay. *Fertil. Steril.* 80(suppl. 3): S41-42.
- Worrlow, K.C. et al. 2004. "Correlation between normal sperm head morphology and sperm binding potential to the sperm HBA. *Fertil. Steril.* 82(Suppl. 2): S95.
- Worrlow, K.C. et al. 2012. Use of hyaluronan in the selection of sperm for intracytoplasmic sperm injection (ICSI); significant improvement in clinical outcomes-multicenter, double blinded and randomized trial. *Hum. Reprod.* 28(2): 306-314.
- Yagci, A. et al. 2010. "Spermatozoa bound to solidstate hyaluronic acid show chromatin structure with high DNA chain integrity: An acridine orange fluorescence study. *J. Androl.* 31(6):566-72.
- Yogev, L. et al. 2010. "Assessing the predictive value of hyaluronan binding ability for the freezability potential of human sperm." *Fertil. Steril.* 93(1):154-8.

## SYMBOLS DEFINITION

Ref #	Symbol	Title	Description	Standard Development Organization
5.1.6		Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.	ISO 15223-1
5.1.5		Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1
5.4.3		Consult instructions for use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1
5.3.4		Keep dry	Indicates a medical device that needs to be protected from moisture.	ISO 15223-1
5.2.7		Non-sterile	Indicates a medical device that has not been subjected to a sterilization process.	ISO 15223-1
5.4.2		Do not re-use	Indicates a medical device that is intended for one use, or use on a single patient during a single procedure.	ISO 15223-1
5.5.1		In vitro diagnostic medical device	Indicates a medical device that is intended to be used as an in vitro diagnostic medical device.	ISO 15223-1
5.2.8		Do not use if package is damaged	Indicates a medical device that should not be used if the package has been broken or damaged.	ISO 15223-1
5.3.7		Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.	ISO 15223-1
5.1.4		Use-by Date	Indicates the date after which the medical device is not to be used.	ISO 15223-1
5.1.2		Authorized representative in the European Community.	Indicates the authorized representative in the European Community	ISO 15223-1
5.1.1		Manufacturer	Indicates the medical device manufacturer, as defined in EU Directive 90/385/EEC, 93/42/EEC and 98/79/EC.	ISO 15223-1
5.5.5		Contains sufficient for <n> tests	Indicates the total number of IVD tests that can be performed with the IVD.	ISO 15223-1
n.a.	<b>Rx Only</b>	Prescription device	Caution: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.	21 CFR 801.109
n.a		Product conforms to the Medical Device Directive 93/42/EEC	Signifies European technical conformity.	n.a.


HBA® is a registered trademark of CooperSurgical, Inc.  
Cell-Vu® is a registered trademarks of CooperSurgical, Inc.


© 2020 CooperSurgical, Inc.

Made in Costa Rica

# CooperSurgical

Phone: (800) 243-2974  
Fax: (800) 262-0105  
**International:**  
Phone: +1 (203) 601-9818  
Fax: +1 (203) 601-4747  
fertility.coopersurgical.com

 CooperSurgical, Inc.  
95 Corporate Drive  
Trumbull, CT 06611 USA

  
CooperSurgical Distribution B.V.  
Celsiusweg 35, 5928PR Venlo,  
Netherlands