English

Quinn's Advantage™
Ca/Mg-Free Medium
with HEPES

For laboratory procedures only; other uses must be qualified by the end user.

Product Description	REF Number	Unit Size
Quinn's Advantage™ Ca/Mg-Free Medium with HEPES	ART-4100	4 x 12 mL

INTENDED USE

This product is intended for in vitro procedures involving the separation of blastomeres from cleaving embryos by disruption of the tight junctional complexes that hold them together.

PRODUCT DESCRIPTION

Ca/Mg-Free Medium with HEPES has been shown in studies (Dumoulin et al, Human Reprod. 1998;13:2880-2883.) to allow easier removal of blastomeres from cleaving embryos by disruption of the tight junctional complexes that hold them together. This medium is a modification of the original Human Tubal Fluid (HTF) medium with HEPES formulation first described by Quinn and colleagues (Fertil Steril. 1984;41:202, 1985;44:493.).

This product contains 10 mg/L of gentamicin, an aminoglycoside antibiotic.

PRECAUTIONS AND WARNINGS

Do not use medium that shows evidence of particulate matter, cloudiness, or is not rose colored.

Single Use:To avoid problems with contamination, handle using aseptic techniques and discard any excess product that remains in the bottle or vial after procedure is completed. Reproductive media products are intended for single use only. Re-use of reproductive media may result in using a product past its labeled expiration date or increase the risk of microbial contamination in a subsequent procedure if the practitioner fails to utilize adequate aseptic techniques.

Use of expired or microbial contaminated product may result in suboptimal conditions to promote fertilization and/or embryo quality during in-vitro culture. These conditions may result in the failure of the embryo to develop properly or to implant, potentially leading to a failed assisted reproductive procedure.

Note: Embryo is considered a general term. More precisely, SAGE™ considers the period of time initiating when a single diploid cell results from the fusion of male and female genome resulting in zygote formation with subsequent development from repeated mitotic divisions forming a solid mass or morula (typically day 4-5) and after which a fluid-filled cavity develops resulting in blastocyst formation (typically day 5-6) ending with embryo implantation that begins the end of the first week and is completed by the end of the second week post conception.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

This product contains the antibiotic gentamicin sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

QUALITY ASSURANCE

Eight-cell MEA tested and passed with 80% or greater blastocyst. USP Endotoxin tested and passed with <1 EU/mL.

A Certificate of Analysis is available for this product.

DIRECTIONS FOR USE

This is the preferred medium for inducing the separation of blastomeres from cleaving embryos done in an atmosphere of air not involving the use of a CO₂ incubator. Cleaving human embryos (4- to 10-cell stage on the third day of development) can be manipulated in this medium for up to 45 minutes with no deleterious effects on subsequent development of the embryo to the blastocyst stage (Dumoulin et al, Human Reprod. 1998;13:2880-2883.). This medium requires the addition of protein or another suitable polymer to stabilize cell membranes during manipulation.

Each laboratory should make its own determination of which medium to use for each particular procedure.

Information on specific aspects of IVF, embryo culture, and cryopreservation is available in our Product Catalog.

STORAGE INSTRUCTIONS AND STABILITY

Store unopened containers at 2-8 °C upon receipt. Warm to incubator (37 °C) temperature prior to use. Do not expose to temperatures greater than 39 °C. The product is stable until the expiration date shown on the label.

- Remove desired volume of product using aseptic procedures.
- Once removed, do not return any volume of product to the original container.
- Do not use if the product becomes discolored, cloudy, turbid, or shows any evidence of microbial contamination.

RELATED PRODUCTS

SAGE In Vitro Fertilization™ has a full line of products for the Reproductive Medicine Specialist. Please call or write for specific information or to receive a copy of our current catalog. For technical questions, or to reach our Customer Service Department, call the SAGE™ Support Line.

Quinn's Advantage™ is a trademark of CooperSurgical, Inc.

**Call the SAGE™ SUPPORT LINE:
In the U.S.: (800) 243-2974 International:
(203) 601-9818**

EXPLANATION OF SYMBOLS

	Catalog Number
	Batch Number
	Use By (year, month, day)
	Do Not Reuse
	Temperature Limitation
	Aseptic Technique Sterilization Membrane Filtered (SAL 10 ⁻³)
	ATTENTION: See instructions for use.
	Authorized Representative in the European Community.
	Product conforms to the Medical Device Directive 93/42/EEC
	Manufacturer

RX ONLY U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

In Vitro Fertilization
a CooperSurgical Company

SAGE™ In Vitro Fertilization, Inc.
a CooperSurgical Company
95 Corporate Drive
Trumbull, CT 06611 USA

ORIGIO a/s
Knardrupvej 2
2760 Måløv
Denmark
www.origio.com

Customer Service:
E-mail: customer.service@origio.com
Tel: +45 46 79 02 02 Fax: +45 46 79 03 02

Français

Quinn's Advantage™
Ca/Mg-Free Medium
with HEPES

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