Instructions for the Use of **global®**

(Catalogue Numbers: LGGG-020, LGGG-050, LGGG-100)

Let the Embryos Choose® "Uninterrupted Culture Medium®"

PRECAUTIONS AND WARNINGS

- Certification: Federal Law (USA) restricts this device to sale by or on the order of a physician (or properly licensed practitioner).
- Caution: The user should read and understand the Instructions for Use, Precautions and Warnings, and be trained in the correct procedure before using **global®** for human embryo culture or transfer.
- Not to be used for injection.
- Do not resterilize.
- Do not use the product if:
 - the product packaging appears damaged or if the seal is broken
 - the expiry date has been exceeded
 - the product becomes discolored, cloudy, or shows evidence of particulate matter
- global®** contains the antibiotic gentamicin sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.
- To avoid problems with contamination, practice aseptic techniques.
- Discard unused medium within 7 days of opening. Do not use after expiry date.

GENERAL INFORMATION

Indications for Use
Culture of human embryos from zygote to blastocyst, embryo transfer.

Storage and Shelf Life

Store at 2-8°C and protected from light. Ten (10) weeks from the date of manufacture.

Composition

A bicarbonate-buffered medium replete with glucose, lactate, pyruvate and all 20 amino acids is optimal to support the growth and development of human embryos in vitro.

Sodium Chloride, Sodium Pyruvate, L-Arginine, L-Threonine, Potassium Chloride, L-Alanine, L-Cysteine, L-Tryptophan, Calcium Chloride, L-Asparagine, L-Histidine, L-Tyrosine, Potassium Phosphate, L-Aspartic Acid, L-Isoleucine, L-Valine, Magnesium Sulfate, L-Glutamic Acid, L-Leucine, Glycyl-L-Glutamine, Sodium Bicarbonate, Glycine, L-Lysine, EDTA, Glucose, L-Proline, L-Methionine, Phenol Red, Sodium Lactate, L-Serine, L-Phenylalanine, Gentamicin Sulfate (*10 µg/ml)

*from therapeutic-grade source material

QUALITY CONTROL SPECIFICATIONS

Assay (performed for each batch)	Specification
Physicochemical Tests	
pH (with 5% CO ₂)	7.2-7.4
Osmolality	260-270 mOsm
Biological Tests	
Endotoxin (LAL)	≤ 0.5 EU/ml
Sterility Test (bacterial and fungal screen, SAL 10 ⁻³)	PASS
Biological Assays	
1-cell Mouse Embryo Assay (% expanded blastocysts at 96 h of culture)	≥ 80%

Special Note on the CO₂ Concentration in the Incubator: In most cases, a 5-7% concentration of CO₂ in the incubator will produce a pH of 7.2 to 7.4 in **global®**. However, the exact concentration of CO₂ required to produce the optimum pH of approximately 7.30 (7.27-7.33) depends on several factors, including the physical characteristics of incubator and the altitude. Consequently, we strongly recommend that each laboratory determine and use the concentration of CO₂ that is required to produce a pH of 7.30 in **global®**.

INSTRUCTIONS FOR USE

The procedures described below have been found to be effective for the preparation of **global®** for culture of human embryos from zygote to blastocyst, and embryo transfer. Every laboratory must define and optimize its own procedures.

After each time the original bottle is opened recap the bottle tightly and store at 2-8°C, protected from light.

Twenty-four (24) hours prior to the use of **global®**, supplement the medium with either Human Serum Albumin (HSA) or LifeGlobal® Protein Supplement to achieve desired (% v/v) of protein supplementation.

- Using a sterile pipette or tip, dispense 25-100 µl droplets or in larger volumes (0.5-1.0 ml) of **global®** supplemented with protein. Cover droplets with appropriate oil.
- Before introducing the embryos, place the culture dishes in the incubator for a minimum of 8 hours to ensure CO₂ and temperature equilibration.

Label each dish with patient information.

- On Day 1, place the zygotes into the equilibrated **global®**. Culture the embryos for 48 h (Day 3, 4-8 cell stage).
- For further culture to blastocyst stage, either a) transfer the embryos to fresh medium under fresh cil and return to the incubator or b) maintain the embryos in the same medium (See Reed et al., 2009; 2010). Note that such uninterrupted culture requires special attention to air quality.

- For transfer on Day 3 (cleavage stage) or Day 5/6 (blastocyst stage) follow general laboratory practice, and transfer to the uterus in 20-30 µl of equilibrated **global®** supplemented with protein.

References
• Reed ML, Hamic A, Thompson DJ and Caperton CL, *Fertil Steril* 92, 1783-6, 2009
• Reed ML, Hamic A, Thompson DJ and Caperton CL, *J. Clin. Embryol.* 13, 33-41, 2010

SYMBOLS

STERILE A	RX Only	REF	LOT	■	■	■
Sterile Usin Asepti Processing Techniques	By Prescription Only	Catalogue Number	Batch Code	Consult Instructions For Use	Manufacturer	Keep Away From Sunlight
2°C	EC REP					
Temperature Limitation	Authorized Representative in the European Community	Use By	GS1 DataMatrix Barcode	Do Not Resterilize	European Conformance (notified body)	Do Not Use If Package is Damaged

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(Références catalogue: LGGG-020, LGGG-050, LGGG-100)

Let the Embryos Choose® "Uninterrupted Culture Medium®"

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