English

**Quinn's Advantage™
Fertilization (HTF) Medium**

For laboratory procedures only; other uses must be qualified by the end user.

Product Description	REF Number	Unit Size
Quinn's Advantage™ Fertilization (HTF) Medium	ART-1020	50 mL
Quinn's Advantage™ Fertilization (HTF) Medium	ART-1021	100 mL

INTENDED USE
This product is intended for in vitro procedures involving the fertilization of human oocytes.

PRODUCT DESCRIPTION
This medium is a modification of Human Tubal Fluid (HTF) that was first described by Quinn and colleagues (Fertil Steril. 1984;41:202, 1985;44:493). Modifications to the original formulation include low phosphate and the addition of citrate, EDTA, selected nonessential amino acids and taurine, alanyl-glutamine as a stable source of glutamine, and the presence of lactate in the form of calcium lactate.

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This product contains 10 mg/L of gentamicin, an aminoglycoside antibiotic.

PRECAUTIONS AND WARNINGS
Do not use medium that shows evidence of particulate matter, cloudiness, or is not rose colored.

Single Use: To avoid problems with contamination, handle using aseptic techniques and discard any excess product that remains in the bottle or vial after procedure is completed.

Reproductive media products are intended for single use only. Re-use of reproductive media may result in using a product past its labeled expiration date or increase the risk of microbial contamination in a subsequent procedure if the practitioner fails to utilize adequate aseptic techniques.

Use of expired or microbial contaminated product may result in suboptimal conditions to promote fertilization and/or embryo quality during in-vitro culture. These conditions may result in the failure of the embryo to develop properly or to implant, potentially leading to a failed assisted reproductive procedure.

The products are aseptically processed and supplied sterile.

Note: Embryo is considered a general term. More precisely, SAGE™ considers the period of time initiating when a single diploid cell results from the fusion of male and female genome resulting in zygote formation with subsequent development from repeated mitotic divisions forming a solid mass or morula (typically day 4-5) and after which a fluid-filled cavity develops resulting in blastocyst formation (typically day 5-6) ending with embryo implantation that begins the end of the first week and is completed by the end of the second week post conception.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

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This product contains the antibiotic gentamicin sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

QUALITY ASSURANCE
One-cell MEA tested and passed with 80% or greater blastocyst. USP Endotoxin tested and passed with <1 EU/mL.

A Certificate of Analysis is available for this product.

DIRECTIONS FOR USE
This is the preferred medium for fertilization procedures in vitro. We currently recommend the use of Quinn's Advantage™ Cleavage Medium (REF #1026/1027) for culture of embryos from Day 1 until Day 3, and Quinn's Advantage™ Blastocyst Medium (REF #1029) for culture from Day 3 to Day 5/6. The media should be supplemented with an exogenous protein source such as Human Serum Albumin (REF #3001/3003) or Serum Protein Substitute (REF #3010/3011).

Each laboratory should make its own determination of which medium to use for each particular procedure.

Information on specific aspects of IVF, embryo culture, and cryopreservation is available in our Product Catalog.

STORAGE INSTRUCTIONS AND STABILITY
Store unopened containers refrigerated at 2 °C to 8 °C. Warm to incubator (37 °C) temperature and equilibrate with desired atmosphere containing 5% CO₂ prior to use. Do not freeze or expose to temperatures greater than 39 °C.

In order to ensure optimal performance, SAGE Media™ strongly recommends measuring the pH of the medium under laboratory working conditions and adjusting the level of CO₂ used to attain the desired pH range for optimal embryo development. The desired pH range for Quinn's Advantage™ Fertilization (HTF) Medium

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is 7.3 ± 0.1. The product is stable until the expiration date shown on the label.

- A. Remove desired volume of product using aseptic procedures.
- B. Once removed, do not return any volume of product to the original container.
- C. Do not use if the product becomes discolored, cloudy, turbid, or shows any evidence of microbial contamination.

RELATED PRODUCTS
ART-1026/1027 Quinn's Advantage™ Cleavage Medium
ART-1029 Quinn's Advantage™ Blastocyst Medium
ART-3001/3003 Human Serum Albumin 100mg/mL in normal saline
ART-3010/3011 Quinn's Advantage™ SPS Serum Protein Substitute
SAGE In Vitro Fertilization™ has a full line of products for the Reproductive Medicine Specialist. Please call or write for specific information or to receive a copy of our current catalog. For technical questions, or to reach our Customer Service Department, call the SAGE™ Support Line.

Quinn's Advantage™ is a trademark of CooperSurgical, Inc.

**Call the SAGE™ SUPPORT LINE:
In the U.S.: (800) 243-2974 International:
(203) 601-9818**

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lečnosti CooperSurgical, Inc.