## English

**Quinn's Advantage™  
Blastocyst Freeze Kit**

For laboratory procedures only; other uses must be qualified by the end user.

| Product Description                      | REF Number | Unit Size |
|------------------------------------------|------------|-----------|
| Quinn's Advantage™ Blastocyst Freeze Kit | ART-8015   | 3 x 12 mL |

**INTENDED USE**

This product is intended for use in freezing blastocysts. This kit is designed to be used in conjunction with the SAGE™ Quinn's Advantage™ Thaw Kit (ART-8016) for warming and recovery of specimens.

**PRODUCT DESCRIPTION**

The components of this kit will allow for the efficient freezing of blastocyst-stage embryos. The components and recommended procedures are the preferred method for improved blastocyst survivability.

This product contains 10 mg/L of genta-micin, an aminoglycoside antibiotic.

**MATERIALS PROVIDED IN THE BLASTOCYST FREEZE KIT**

- 1 x 12 mL vial of 5% Glycerol Freezing Medium (REF #8009-12) with 12 mg/mL Human Serum Albumin
- 1 x 12 mL vial of 9% Glycerol with 0.2 M Sucrose Freezing Medium (REF #8011-12) with 12 mg/mL Human Serum Albumin
- 1 x 12 mL vial of Freeze/Thaw Diluent Solution (REF #8013-12) with 12 mg/mL Human Serum Albumin

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**PRECAUTIONS AND WARNINGS**

Do not use medium that shows evidence of particulate matter, cloudiness, or is not rose colored.

The 5% Glycerol Freezing Medium (ART-8009-12) component in this kit contains 12 mg/mL human serum albumin.

The 9% Glycerol with 0.2 M Sucrose Freezing Medium (ART-8011-12) component in this kit contains 12 mg/mL human serum albumin.

The Freeze/Thaw Diluent Solution (ART-8013-12) component in this kit contains 12 mg/mL human serum albumin.

**Caution:** All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when testing for antibodies to HIV-1/HIV-2, HCV and non-reactive for HBsAg, HCV RNA and HIV-1 RNA. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents. Donors of the source material have been screened for Creutzfeldt-Jacobs disease (CJD). Based on effective donor screening and product manufacturing processes, it carries an extremely remote risk for transmission of viral diseases. A theoretical risk for transmission of CJD is also considered extremely remote. No cases of transmission of viral diseases or CJD have ever been identified for albumin.

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There

are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes.

**Single use:** To avoid problems with contamination, handle using aseptic techniques and discard any excess product that remains in the bottle or vial after procedure is completed.

Reproductive media products are intended for single use only. Re-use of reproductive media may result in using a product past its labeled expiration date or increase the risk of microbial contamination in a subsequent procedure if the practitioner fails to utilize adequate aseptic techniques.

Use of expired or microbial contaminated product may result in suboptimal conditions to promote fertilization and/or embryo quality during in-vitro culture. These conditions may result in the failure of the embryo to develop properly or to implant, potentially leading to a failed assisted reproductive procedure.

**Note:** Embryo is considered a general term. More precisely, SAGE™ considers the period of time initiating when a single diploid cell results from the fusion of male and female genome resulting in zygote formation with subsequent development from repeated mitotic divisions forming a solid mass or morula (typically day 4-5) and after which a fluid-filled cavity develops resulting in blastocyst formation (typically day 5-6) ending with embryo implantation that begins the end of the first week and is completed by the end of the second week post conception.

**Caution:** U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

This product contains the antibiotic genta-micin sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

\*Vials, eg, 1.2 mL plastic cryovials, also can be

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**QUALITY ASSURANCE**

One-cell MEA tested and passed with 80% or greater blastocyst. USP Endotoxin tested and passed with <1 EU/mL.

A Certificate of Analysis is available for this product.

**DIRECTIONS FOR USE FOR****CRYOPRESERVATION OF BLASTOCYSTS**

Controlled hyperstimulation of women undergoing IVF or GIFT produces, on average, 10 to 12 mature oocytes for insemination. It is prudent to replace only a limited number of the resulting embryos, as multiple pregnancies can arise if too many embryos are replaced. Therefore, the majority of patients will have supernumerary embryos. These embryos can be cryopreserved and stored for later use, thus avoiding the necessity of the couple to undertake another stimulated cycle to recover more oocytes for IVF.

The major cause of cell damage during cryopreservation is the formation of intracellular ice during freezing and thawing. By using cryopreservatives, controlling the rates of freezing and thawing, and carefully diluting the cryopreservative from the embryo after thawing, methods have been developed that allow 80% or more of frozen-thawed embryos to survive and be replaced into the reproductive tract of the woman who produced the oocytes or a genetically unrelated recipient. Place fully expanded blastocysts into 1.0 mL of Freeze/Thaw Diluent Solution (REF #8013-12) at 37 °C for 5 minutes. Transfer to 1.0 mL of the 5% Glycerol Freezing Medium (REF #8009-12) at 37 °C for 10 minutes. Transfer to 1.0 mL of the 9% Glycerol + 0.2 M Sucrose Freezing Medium (REF #8011-12) at 37 °C and load into straws\* containing the same solution; hold at 37 °C in the 9% Glycerol + 0.2 M Sucrose Freezing Medium for a total of 10 minutes before cooling is initiated.

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used and have been found to be preferable by some laboratories.

It is important to make sure that the embryos are well mixed with the cryoprotectant solutions. This can be accomplished by pipetting the embryos up and down in the solution several times after adding them to the cryoprotectant solution. It is also recommended that the media be covered with Sterile Oil for Tissue Culture (REF #4008) during use to minimize evaporation of water and a subsequent change in osmolality of the solutions.

**COOLING PROTOCOL**

Embryos are usually frozen in straws, but 1.2 mL vials can be used as an alternative. The embryos are taken from a starting temperature of 37 °C to -6 °C at 2 °C/min. They are then seeded manually and held at -6 °C for a total of 10 to 15 minutes before being cooled at about 0.3 °C/min to around -35 °C. They are then transferred to a storage tank of liquid nitrogen.

**Each laboratory should make its own determination of which medium to use for each particular procedure.**

Information on specific aspects of IVF, embryo culture, and cryopreservation is available in our Product Catalog.

**STORAGE INSTRUCTIONS AND STABILITY**

Store unopened containers refrigerated at 2 °C to 8 °C. Bring to 37 °C prior to use. Do not freeze or expose to temperatures greater than 39 °C. The product is stable until the expiration date shown on the label.

A. Remove desired volume of product using aseptic procedures.

B. Once removed, do not return any volume of product to the original container.

C. Do not use if the product becomes discolored, cloudy, turbid, or shows any evidence of microbial contamination.

**RELATED PRODUCTS**

ART-8016 Quinn's Advantage™ Thaw Kit  
ART-4008P Oil for Tissue Culture

SAGE In Vitro Fertilization™ has a full line of products for the Reproductive Medicine Specialist. Please call or write for specific information or to receive a copy of our current catalog. For technical questions, or to reach our Customer Service Department, call the SAGE™ Support Line.

Call the SAGE™ SUPPORT LINE: In the U.S.: (800) 243-2974 International: (203) 601-9818

**EXPLANATION OF SYMBOLS**

|  |                                                                                                                       |
|--|-----------------------------------------------------------------------------------------------------------------------|
|  | Catalog Number                                                                                                        |
|  | Batch Number                                                                                                          |
|  | Use By (year, month, day)                                                                                             |
|  | Do Not Reuse                                                                                                          |
|  | Temperature Limitation                                                                                                |
|  | Aseptic Technique Sterilization Membrane Filtered (SAL 10 <sup>-3</sup> )                                             |
|  | <b>ATTENTION:</b> See instructions for use.                                                                           |
|  | Authorized Representative in the European Community.                                                                  |
|  | Product conforms to the Medical Device Directive 93/42/EEC                                                            |
|  | Manufacturer                                                                                                          |
|  | U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner). |

**SAGE** In Vitro Fertilization  
*a CooperSurgical Company*

SAGE™ In Vitro Fertilization, Inc.  
a CooperSurgical Company  
95 Corporate Drive  
Trumbull, CT 06611 USA

**Customer Service:**  
E-mail: customer.service@origio.com  
Tel: +45 46 79 02 02 Fax: +45 46 79 03 02

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