MediCult IVM® System

Product No.:

8221

CE

2797

Languages:

bg

cs

da

de

el

en

es

et

fr

hr

hu

is

it

kk

lt

lv

mk

nl

no

pt

ro

ru

sk

sl

sv

tr

uk

Customer Service:
E-mail: customer.service@origio.com
Tel: +45 46 79 02 02 • Fax: +45 46 79 03 02

ORIGIO a/s
Knardrupvej 2, DK-2760 Måløv, Denmark
www.origio.com
Tel.: +45 46 79 02 00 Fax: +45 46 79 03 00

a CooperSurgical Company

Ingredientia
8201 LAG Medium – Vial 1
Albumini humani solutio; Aqua; CaCl₂;
Gentamicinum; Glucosum; KCl; MgSO₄;
NaCl; NaH₂PO₄; NaHCO₃; Pyruvas;
Phenolsulfonphthaleinum; SSR

8220 IVM® Medium – Vial 2
Acetas; Acida aminica: Ala, Arg, Asp, Cys,
Gln, Glu, Gly, His, Ile, Leu, Lys, Met, Phe,
Pro, Ser, Thr, Trp, Tyr, Val, Cystinum,
Glutathionum, Hydroxyprolinum; Adeninium;
Adenosini phosphas; Adenosini triphosphas;
Aqua; CaCl₂; Cholesterolum; Desoxy-D-
ribosum; D-ribosum; Fe(NO₃)₃; Gentamicinum;
Glucosum; Guaninium; Hypoxanthinum;
KCl; MgSO₄; NaCl; NaH₂PO₄; NaHCO₃;
Phenolsulfonphthaleinum; Polysorbatum 80;
Pyruvas; Thyminium; Uracilium; Vitamini:
A, B₁, B₂, B₃, B₅, B₆, C, D₂, E, Acidum
4-aminobenzoicum, Acidum nicotinicum,
Cholinum, Inositolum, Menadionum,
Nicotinamidum, Pyridoxalum, Pyridoxinum;
Xanthinum

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en - english

For pre-incubation and maturing of immature oocytes.

This product is for ART treatment, whether the cause of infertility is male or female. The product should only be used by professionals trained in ART treatment.

Package

8201 LAG Medium (Vial 1) 8220 IVM® Medium (Vial 2)

Pack size

1 x Vial 1 (10 ml) 3 x Vial 2 (10 ml)

Contains

Human albumin solution (HAS) (only Vial 1) Recombinant human insulin (only Vial 1) Gentamicin sulphate 10 µg/ml

Quality control testing

Sterility tested (Ph.Eur., USP) Osmolality tested (Ph.Eur., USP) pH tested (Ph.Eur., USP) Endotoxin tested ≤ 0.1 EU/ml (Ph.Eur., USP) HSA analysis (Ph.Eur., USP) Mouse Embryo Assay (MEA) tested
Note: The results of each batch are stated on a Certificate of Analysis, which is available at www.origio.com.

Storage instructions and stability

The products are aseptically processed and supplied sterile. Store in original container at 2-8°C, protected from light. Do not freeze. Discard excess (unused) media following warming. IVM® Medium is provided in vials intended for single use. LAG Medium is to be used within 7 days after opening. When stored as directed by the manufacturer the product is stable until the expiry date shown on the label.

Precautions and warnings

Do not use the product if:

- Product packaging appears damaged or if the seal is broken.
- Expiry date has been exceeded.
- The product becomes discoloured, cloudy, turbid, or shows any evidence of microbial contamination.

Caution: The product is for single use and is not to be reused due to risk of contamination.

Caution: All blood products should be treated as potentially infectious. The source material used to manufacture this product was tested and found non-reactive for HBsAg and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.

Note: Please note the need for traceability of this product. In addition, national legal requirements for this field may exist in your

country.

Note: Only to be used in combination with other devices intended for the particular purpose.

Note: Dispose of the device in accordance with local regulations for disposal of medical devices.

Instructions for use

- Pre-equilibrate 3 ml LAG Medium (Vial 1) and 10 ml of IVM® Medium (Vial 2) in a 5-6% CO₂ environment at 37°C for a minimum of 12 hours.
- After oocyte pick up, store the immature oocytes in LAG Medium in a 5-6% CO₂ environment at 37°C for 2-3 hours prior to transfer to the final IVM® Maturation Medium.
- Preparation of final IVM® Maturation Medium: 9 ml IVM® Medium (Vial 2) 1 ml patient's own serum 10 µl human chorionic gonadotropin (hCG) solution (100 mIU/ml) 100 µl human follicle stimulating hormone (FSH) solution (75 mIU/ml).
- Transfer the oocytes to the final IVM® Maturation Medium and incubate in a 5-6% CO₂ environment at 37°C for 28-32 hours.

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en - english	
	
Caution: This product is intended for assisted fertilization, regardless of the cause of infertility. The product should only be used by trained personnel.	
Packaging 8201 LAG Medium (Vial 1) 8220 IVM® Medium (Vial 2)	
Packaging size 1 x Vial 1 (10 ml) 3 x Vial 2 (10 ml)	
Contents Human albumin solution (HAS) (Vial 1 only) Recombinant human insulin (Vial 1 only) Gentamicin sulfate 10 µg/ml	
Quality control testing Sterility test (Ph.Eur., USP) Osmolality test (Ph.Eur., USP) pH test (Ph.Eur., USP) Test for endotoxin ≤ 0.1 EU/ml (Ph.Eur., USP) HSA analysis (Ph.Eur., USP) Embryo analysis (MEA), tested	
Notes: Results of each batch are given on an analysis certificate available at www.origo.com .	

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