




WALLACE® EMBRYO TRANSFER

Wallace Classic


Embryo Transfer Catheters



**Hand crafted soft
inner catheter**



**Smooth tip transition
with outer sheath**



**Inner catheter
position in uterus**

Wallace Classic

Since the very first IVF case, the Wallace Classic range continues to offer a soft-profiled tip

- Soft material and smooth tip minimizes trauma and maximizes success rates
- Formable outer sheath to suit different patient anatomies
- Centimeter markings to show outer sheath depth and inner lumen protrusion for confidence of placement within the uterus
- Trial transfer catheters allow for confident positioning of the outer sheath before introducing the embryo loaded inner. Clinicians have found this technique to have a positive impact on pregnancy rates¹
- The Classic range is fully compatible with the SureView range

Product Configurations

Classic Trial Transfer Catheter (18cm/23cm)
Classic Embryo Transfer Catheter (18cm/23cm)
Intermediate Embryo Transfer Catheter (18cm/23cm)
Malleable Stylet - for difficult transfer (18cm/23cm)

1. Neithardt, A.B et al (2005). Embryo afterloading: a refinement in embryo transfer technique that may increase clinical pregnancy. Fertility and Sterility. Vol. 83, No.3.

Technical Information

Wallace Classic / Wallace SureView	
Inner catheter (ID)	0.76mm
Inner catheter (OD)	1.52mm
Outer sheath (OD)	2.3mm
Volume of catheter not including hub (18cm catheter)	0.081ml (81µl)
Volume of catheter not including hub (23cm catheter)	0.104ml (104µl)
Volume per cm length of catheter	0.0045ml (4.5µl)
Suggested loading volume	0.5 to 1cm = 2.25 to 4.5µl
Stylet (OD)	1.52mm
Stylet protrusion from outer sheath	2.0mm

Product Compatibility and Order Codes

Compatible with:	Wallace Classic			
	1816 Qty: 10 Soft 18cm Classic	1816N Qty: 10 Soft 23cm Classic	1816S Qty: 10 Intermediate 18cm Classic	2316 Qty: 10 Intermediate 23cm Classic
1816ST Quantity: 10 Rigid stylet for introducing the 18cm	•		•	
1816NST Quantity: 10 Rigid stylet for introducing the 23cm		•		•
TT1816 Quantity: 10 Trial transfer catheter blind end (18cm) Classic catheter	•		•	
TT1816N Quantity: 10 Trial transfer catheter blind end (23cm) Classic catheter		•		•

• = Compatible

CE marked in compliance with EU medical device directive