




EMBRYO TRANSFER BY WALLACE®

Wallace Classic


Embryo Transfer Catheters



Hand crafted soft inner catheter



Smooth tip transition with outer sheath



Inner catheter position in uterus



CooperSurgical Fertility Companies

ORIGIO a/s | origio.com

CooperGenomicsSM
a CooperSurgical company

globalsupport@coopergenomics.com
coopergenomics.com

Since the very first IVF case, the Wallace Classic range continues to offer a soft-profiled tip

- Soft material and smooth tip minimizes trauma and maximizes success rates
- Formable outer sheath to suit different patient anatomies
- Centimeter markings to show outer sheath depth and inner lumen protrusion for confidence of placement within the uterus
- Trial transfer catheters allow for confident positioning of the outer sheath before introducing the embryo loaded inner. Clinicians have found this technique to have a positive impact on pregnancy rates¹
- The Classic range is fully compatible with the SureView range

Product Configurations
Classic Trial Transfer Catheter (18cm/23cm)
Classic Embryo Transfer Catheter (18cm/23cm)
Intermediate Embryo Transfer Catheter (18cm/23cm)
Malleable Stylet - for difficult transfer (18cm/23cm)

1. Neihardt, A.B et al (2005). Embryo afterloading: a refinement in embryo transfer technique that may increase clinical pregnancy. Fertility and Sterility, Vol. 83, No.3.

Technical Information	Wallace Classic / Wallace SureView
Inner catheter (ID)	0.76mm
Inner catheter (OD)	1.52mm
Outer sheath (OD)	2.3mm
Volume of catheter not including hub (18cm catheter)	0.081ml (81µl)
Volume of catheter not including hub (23cm catheter)	0.104ml (104µl)
Volume per cm length of catheter	0.0045ml (4.5µl)
Suggested loading volume	0.5 to 1cm = 2.25 to 4.5µl
Stylet (OD)	1.52mm
Stylet protrusion from outer sheath	2.0mm

Product Compatibility and Order Codes:

	Wallace Classic			
	1816 Quantity: 10 Soft 18cm Classic	1816N Quantity: 10 Soft 23cm Classic	1816S Quantity: 10 Intermediate 18cm Classic	2316 Quantity: 10 Intermediate 23cm Classic
Compatible with:				
1816ST Quantity: 10 Rigid stylet for introducing the 18cm	✓		✓	
1816NST Quantity: 10 Rigid stylet for introducing the 23cm		✓		✓
TT1816 Quantity: 10 Trial transfer catheter blind end (18cm) Classic catheter	✓		✓	
TT1816N Quantity: 10 Trial transfer catheter blind end (23cm) Classic catheter		✓		✓

✓ = Compatible

CE marked in compliance with EU medical device directive



CooperSurgical Fertility Companies

ORIGIO a/s | origio.com

CooperGenomicsSM
a CooperSurgical company

globalsupport@coopergenomics.com
coopergenomics.com